

Olerup SSP® HLA-B*57:01

Product number: 101.572-12 – including *Taq* polymerase
101.572-12u - without *Taq* polymerase
Lot number: 12K
Expiry date: 2012-June-01
Number of tests: 12
Number of wells per test: 11

CHANGES COMPARED TO THE PREVIOUS HLA-B*57:01 LOT (44G):

Well	5'-primer	3'-primer	rationale
3	-	Added	Primer added for the B*57:27 allele.
4	Added	Added	Primer pairs added for the B*57:20 and 57:29 alleles.
5	Added	Added	Primers added for the B*57:23 and 57:26 alleles.
9	New	New	New primer pair for the B*57:21 allele.
10	New	New	New primer pair for the B*57:13 and 57:22 alleles.
11	New	New	Negative control.

THE NUMBER OF WELLS has been increased from 8 to 11.

ALLELE COVERAGE:

B*57:01:01 to B*57:01:07, i.e. all the currently recognized HLA-B*57:01 alleles, gives rise to a unique amplification pattern; www.ebi.ac.uk/imgt/hla, 2010-April-01, release 3.0.0.

RESOLUTION IN HLA-B*57:01 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-B*57:01 SUBTYPINGS BY NON-HLA-B*57 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.