

Olerup SSP® HLA-B*57:01

Product number: 101.572-12 – including *Taq* polymerase
101.572-12u - without *Taq* polymerase
Lot number: 08M
Expiry date: 2013-October-01
Number of tests: 12
Number of wells per test: 12+1

CHANGES COMPARED TO THE PREVIOUS HLA-B*57:01 LOT (12K):

Well	5'-primer	3'-primer	rationale
5	Modified, added	-	Improved specificity, primer added for the B*57:34 allele.
9	Added	-	Primers added for the B*57:33 and 57:40 alleles.
10	Added	Added	Primer pairs added for the B*57:37, 57:41 and 57:43 alleles.
11	New, removed	New, removed	New primer pair for improved allelic resolution, negative control moved to well 13.
12	New	New	New primer pair for the B*57:35, 57:36 and 57:38 alleles.
13	Moved	Moved	Negative control moved from well 11.

THE NUMBER OF WELLS has been increased from 11 to 13.

ALLELE COVERAGE:

B*57:01:01 to B*57:01:10, i.e. all the currently recognized HLA-B*57:01 alleles, gives rise to a unique amplification pattern; www.ebi.ac.uk/imgt/hla, 2011-January-14, release 3.3.0.

RESOLUTION IN HLA-B*57:01 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-B*57:01 SUBTYPINGS BY NON-HLA-B*57 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.