

101.572-12 – including *Taq* polymerase101.572-12u – without *Taq* polymeraseLot No.: **86M****Olerup SSP® HLA-B\*57:01**

**Product number:** 101.572-12 – including *Taq* polymerase  
101.572-12u - without *Taq* polymerase  
**Lot number:** 86M  
**Expiry date:** 2014-May-01  
**Number of tests:** 12  
**Number of wells per test:** 12+1

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*57:01 LOT (08M):**

Well	5'-primer	3'-primer	rationale
8	Added	-	Primer added for the B*57:44 allele.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

B\*57:01:01 to B\*57:01:11, i.e. all the currently recognized HLA-B\*57:01 alleles, gives rise to a unique amplification pattern; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2011-April-08, release 3.4.0.

**RESOLUTION IN HLA-B\*57:01 HOMO- AND HETEROZYGOTES:**

Good.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*57:01 SUBTYPINGS BY NON-HLA-B\*57 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.