

HLA-B*57:01
101.572-12 – including *Taq* polymerase
101.572-12u – without *Taq* polymerase
Lot No.: 8L3

Olerup SSP[®] HLA-B*57:01

Product number: 101.572-12 – including *Taq* polymerase
 101.572-12u – without *Taq* polymerase
Lot number: 8L3
Expiry date: 2025-01-01
Number of tests: 12
Number of wells per test: 23+1

CHANGES COMPARED TO THE PREVIOUS HLA-B*57:01 LOT (4K0):

Well	5'-primer	3'-primer	rationale
22	Modified	-	Modified 5'-primer for improved HLA-specific amplification.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*57:01:01 to B*57:01:43, i.e. all the currently recognized HLA-B*57:01 alleles, gives rise to a unique amplification pattern¹; www.ebi.ac.uk/imgt/hla, 2020-October-15, release 3.42.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HLA-B*57:01 HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-B*57:01 SUBTYPINGS BY NON-HLA-B*57 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.