

HLA-B\*58  
101.568-06 – licensed for PCR  
101.568-06u – not licensed for PCR  
Lot No.: Y35

1

## Release Note - Olerup SSP™ HLA-B\*58 – Lot No. Y35

Product number: 101.568-06 – licensed for PCR  
101.568-06u - not licensed for PCR  
Lot number: Y35  
Expiry date: 2009-September-01  
Number of tests: 6  
Number of tubes per test: 16

### CHANGES COMPARED TO THE PREVIOUS HLA-B\*58 LOT (V83):

Tube	5'-primer	3'-primer	rationale
4	-	Added	Primer added for the B*5815 allele.
9	Added	Added	Primer pair added for the B*5814 allele.

THE NUMBER OF WELLS is unchanged.

### ALLELE COVERAGE:

B\*5801 to B\*5815, i.e. all the currently recognized HLA-B\*58 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2007-October-05, release 2.19.0.

### RESOLUTION IN HLA-B\*58 HOMO- AND HETEROZYGOTES:

Good.

### INFLUENCE ON THE INTERPRETATION OF HLA-B\*58 SUBTYPINGS BY NON-HLA-B\*58 ALLELES:

None frequently occurring.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.