

Lot No.: **37G**

Olerup SSP[®] HLA-B*59

Product number: 101.554-06 – including *Taq* polymerase
101.554-06u – without *Taq* polymerase
Lot number: 37G
Expiry date: 2011-October-01
Number of tests: 6
Number of wells per test: 8

CHANGES COMPARED TO THE PREVIOUS HLA-B*59 LOT 02F):

Well	5'-primer	3'-primer	rationale
7	New	New	New primer pair for the B*5904 allele.
8	New	New	New primer pair for the B*5905 allele.

THE NUMBER OF WELLS has been increased from 6 to 8.

ALLELE COVERAGE:

B*5901 to B*5905, i.e. all the currently recognized HLA-B*59 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-July-17, release 2.26.0.

RESOLUTION IN HLA-B*59 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*59 SUBTYPINGS BY NON-HLA-B*59 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.