

Olerup SSP[®] HLA-B*59

Product number: 101.554-06 – including *Taq* polymerase
101.554-06u – without *Taq* polymerase
Lot number: 95K
Expiry date: 2013-August-01
Number of tests: 6
Number of wells per test: 8

CHANGES COMPARED TO THE PREVIOUS HLA-B*59 LOT 37G):

The HLA-B*59 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*59:01 to B*59:05, i.e. all the currently recognized HLA-B*59 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.

RESOLUTION IN HLA-B*59 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*59 SUBTYPINGS BY NON-HLA-B*59 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.