

101.550-06 – including *Taq* polymerase101.550-06u – without *Taq* polymeraseLot No.: **45F****Olerup SSP<sup>®</sup> HLA-B\*67**

**Product number:** 101.550-06 – including *Taq* polymerase  
101.550-06u – without *Taq* polymerase  
**Lot number:** 45F  
**Expiry date:** 2011-January-01  
**Number of tests:** 6  
**Number of wells per test:** 4

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*67 LOT (Y34):**

Well	5'-primer	3'-primer	rationale
4	-	Modified	Increased specificity of specific primer pair.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

B\*6701 to B\*6702, i.e. all the currently recognized HLA-B\*67 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-October-10, release 2.23.0.

**RESOLUTION IN HLA-B\*67 HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*67 SUBTYPINGS BY NON-HLA-B\*67 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.