

Olerup SSP[®] HLA-B*67

Product number: 101.550-06 – including *Taq* polymerase
101.550-06u – without *Taq* polymerase
Lot number: 16K
Expiry date: 2012-May-01
Number of tests: 6
Number of wells per test: 4

CHANGES COMPARED TO THE PREVIOUS HLA-B*67 LOT (45F):

The HLA-B*67 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*67:01 to B*67:02, i.e. all the currently recognized HLA-B*67 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-April-01, release 3.0.0.

RESOLUTION IN HLA-B*67 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*67 SUBTYPINGS BY NON-HLA-B*67 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.