

101.550-06 – including *Taq* polymerase101.550-06u – without *Taq* polymeraseLot No.: **49M****Olerup SSP® HLA-B\*67**

**Product number:** 101.550-06 – including *Taq* polymerase  
101.550-06u – without *Taq* polymerase  
**Lot number:** 49M  
**Expiry date:** 2014-January-01  
**Number of tests:** 6  
**Number of wells per test:** 6

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*67 LOT (16K):**

Well	5'-primer	3'-primer	rationale
1	Modified	Modified	Improved specificity of primer pair.
5	New	New	Improved allelic resolution.
6	New	New	New primer pair for the B*67:03 allele.

**THE NUMBER OF WELLS** is increased from 4 to 6.

**ALLELE COVERAGE:**

B\*67:01 to B\*67:03, i.e. all the currently recognized HLA-B\*67 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2011-April-08, release 3.4.0.

**RESOLUTION IN HLA-B\*67 HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*67 SUBTYPINGS BY NON-HLA-B\*67 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.