

Olerup SSP[®] HLA-B*67

Product number:	101.550-06 – including <i>Taq</i> polymerase 101.550-06u – without <i>Taq</i> polymerase
Lot number:	76R
Expiry date:	2015-May-01
Number of tests:	6
Number of wells per test:	6

CHANGES COMPARED TO THE PREVIOUS HLA-B*67 LOT (49M):

The HLA-B*67 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*67:01 to B*67:03, i.e. all the currently recognized HLA-B*67 alleles, will be amplified by the primers in the HLA-B*67 kit; www.ebi.ac.uk/imgt/hla, 2012-October-17, release 3.10.0.

The HLA- B*67 kit enables separation of the confirmed HLA-B*67 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA- B*67 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

The HLA-B*67 subtyping kit cannot distinguish the silent mutations in the B*67:01:01 and B*67:01:02 alleles.

RESOLUTION IN HLA-B*67 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*67 SUBTYPINGS BY NON-HLA-B*67 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.