

HLA-B*67
101.550-06 – licensed for PCR
101.550-06u – not licensed for PCR
Lot No.: Y34

1

Release Note - *Olerup SSP*TM HLA-B*67 – Lot No. Y34

Product number:	101.550-06 – licensed for PCR 101.550-06u – <u>not</u> licensed for PCR
Lot number:	Y34
Expiry date:	2009-September-01
Number of tests:	6
Number of tubes per test:	4

CHANGES COMPARED TO THE PREVIOUS HLA-B*67 LOT (X02):

The HLA-B*67 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*6701 to B*6702, i.e. all the currently recognized HLA-B*67 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-October-05, release 2.19.0.

RESOLUTION IN HLA-B*67 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*67 SUBTYPINGS BY NON-HLA-B*67 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.