

## **Olerup SSP<sup>®</sup> HLA-B\*07 Add-on**

|                           |   |
|---------------------------|---|
| Product number:           | 101.852-12 – including <i>Taq</i> polymerase<br>101.852-12u – without <i>Taq</i> polymerase |
| Lot number:               | 23N   |
| Expiry date:              | 2014-August-01  |
| Number of tests:          | 12  |
| Number of wells per test: | 5   |

THE NUMBER OF WELLS is 5 in an 8 well plate.

### **ALLELE COVERAGE:**

The HLA-B\*07:02, 07:05, 07:06 and 07:61 alleles give different patterns in the B\*07 Add-on subtyping kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2012-January-12, release 3.7.0.

<sup>1</sup>The HLA-B\*07 Add-on kit cannot distinguish these silent mutations: the B\*07:02:01-07:02:29 and the B\*07:05:01-07:05:06 alleles.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF B\*07 ADD-ON SUBTYPINGS BY NON- HLA-B\*07 ALLELES:**

None of importance.