

Olerup SSP[®] HLA-B*07 Add-on

Product number:	101.852-12 – including <i>Taq</i> polymerase 101.852-12u – without <i>Taq</i> polymerase
Lot number:	7L6
Expiry date:	2024-11-01
Number of tests:	12
Number of wells per test:	6+1

CHANGES COMPARED TO THE PREVIOUS HLA-B*07 ADD-ON LOT (8G4):

The HLA-B*07 Add-on primer set is unchanged compared to the previous *Olerup SSP[®] HLA-B*07 Add-on (Lot No. 8G4)*.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The HLA-B*07:02, 07:05, 07:06, 07:61 and 07:161N alleles give different patterns in the B*07 Add-on subtyping kit¹; www.ebi.ac.uk/imgt/hla, 2020-July-13, release 3.41.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF B*07 ADD-ON SUBTYPINGS BY NON- HLA-B*07 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.