

Olerup SSP® HLA-B*73

Product number:	101.555-06 – including <i>Taq</i> polymerase 101.555-06u – without <i>Taq</i> polymerase
Lot number:	52K
Expiry date:	2013-April-01
Number of tests:	6
Number of wells per test:	3

THE NUMBER OF WELLS is 3 wells in an 8-well plate.

ALLELE COVERAGE:

B*73:01 to B*73:02, i.e. all the currently recognized HLA-B*73 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-July-01, release 3.1.0.

RESOLUTION IN HLA-B*73 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*73 SUBTYPINGS BY NON-HLA-B*73

ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.