

Olerup SSP® HLA-B*73

Product number:	101.555-06 – including <i>Taq</i> polymerase 101.555-06u – without <i>Taq</i> polymerase
Lot number:	11R
Expiry date:	2015-February-01
Number of tests:	6
Number of wells per test:	3

CHANGES COMPARED TO THE PREVIOUS HLA-B*73 LOT (52K):

The HLA-B*73 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*73:01 to B*73:02, i.e. all the currently recognized HLA-B*73 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2012-July-12, release 3.9.0.

The HLA-B*73 kit enables separation of the confirmed HLA-B*73 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-B*73 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles

RESOLUTION IN HLA-B*73 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*73 SUBTYPINGS BY NON-HLA-B*73 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.