

101.555-06 – including *Taq* polymerase101.555-06u – without *Taq* polymeraseLot No.: **2F1****Olerup SSP<sup>®</sup> HLA-B\*73**

**Product number:** 101.555-06 – including *Taq* polymerase  
101.555-06u – without *Taq* polymerase

**Lot number:** 2F1

**Expiry date:** 2019-10-01

**Number of tests:** 6

**Number of wells per test:** 3+1

The HLA-B\*73 primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

B\*73:01 to B\*73:02, i.e. all the currently recognized HLA-B\*73 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2016-October-14, release 3.26.0<sup>1</sup>.

The HLA-B\*73 kit enables separation of the confirmed HLA-B\*73 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-B\*73 kit also enables identification of polymorphisms in exons outside of the region encoding the peptide binding domain and of null and alternatively expressed alleles.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

**RESOLUTION IN HLA-B\*73 HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*73 SUBTYPINGS BY NON-HLA-B\*73****ALLELES:**

None.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.