

101.555-06 – including *Taq* polymerase101.555-06u – without *Taq* polymeraseLot No.: **2H4****Olerup SSP[®] HLA-B*73**

Product number: 101.555-06 – including *Taq* polymerase
101.555-06u – without *Taq* polymerase

Lot number: 2H4

Expiry date: 2023-03-01

Number of tests: 6

Number of wells per test: 3+1

CHANGES COMPARED TO THE PREVIOUS HLA-B*73 LOT (2F1):

The HLA-B*73 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*73:01 to B*73:02, i.e. all the currently recognized HLA-B*73 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2019-January-23, release 3.35.0¹.

The HLA-B*73 kit enables separation of the confirmed HLA-B*73 alleles as listed in the IMGT/HLA database. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-B*73 kit also enables identification of null and alternatively expressed alleles.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HLA-B*73 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*73 SUBTYPINGS BY NON-HLA-B*73**ALLELES:**

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to shelf-life extension.