

HLA-B*73
101.555-06 – including *Taq* polymerase
101.555-06u – without *Taq* polymerase
Lot No.: 2R1

Olerup SSP® HLA-B*73

Product number:	101.555-06 – including <i>Taq</i> polymerase 101.555-06u – without <i>Taq</i> polymerase
Lot number:	2R1
Expiry date:	2026-07-01
Number of tests:	6
Number of wells per test:	3+1

CHANGES COMPARED TO THE PREVIOUS HLA-B*73 LOT (2H4):

The HLA-B*73 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*73:01 to B*73:03, i.e. all the currently recognized HLA-B*73 alleles, will be amplified by the primers in the HLA-B*73 subtyping kit; www.ebi.ac.uk/imgt/hla, 2022-April-19, release 3.48.0¹.

The HLA-B*73 kit enables separation of the confirmed HLA-B*73 alleles as listed in the IMGT/HLA database 3.35.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-B*73 kit also enables identification of many null and alternatively expressed alleles.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HLA-B*73 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*73 SUBTYPINGS BY NON-HLA-B*73 ALLELES:

None.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.



0197