

Olerup SSP[®] HLA-B*78

Product number: 101.551-06 – including *Taq* pol.
101.551-01u – without *Taq* pol.
Lot number: 66F
Expiry date: 2011-March-01
Number of tests: 6
Number of tubes per test: 13

CHANGES COMPARED TO THE PREVIOUS HLA-B*78 LOT (Y33):

Well	5'-primer	3'-primer	rationale
5	-	-	Exchanged positive control primer pair.
10	Modified	-	Modified 5'-primer for increased yield of specific PCR product.
13	New	New	New primer pair for the B*7806 allele.

THE NUMBER OF WELLS has been increased from 12 to 13.

ALLELE COVERAGE:

B*7801 to B*7806, i.e. all the currently recognized HLA-B*78 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2009-January-16, release 2.24.0.

RESOLUTION IN HLA-B*78 HOMO- AND HETEROZYGOTES:

The B*7801,7801, B*7801,7802 and B*7802,7802 genotypes give rise to unique amplification patterns.

INFLUENCE ON THE INTERPRETATION OF HLA-B*78 SUBTYPINGS BY NON-HLA-B*78 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.