

Olerup SSP[®] HLA-B*78

Product number: 101.551-06 – including *Taq* pol.
101.551-06u – without *Taq* pol.
Lot number: 48K
Expiry date: 2012-September-01
Number of tests: 6
Number of tubes per test: 15

CHANGES COMPARED TO THE PREVIOUS HLA-B*78 LOT (66F):

Well	5'-primer	3'-primer	rationale
10	-	Added	Increased yield of specific PCR product.
14	New	New	Improved allelic resolution.
15	New	New	New primer pair for the B*78:07 allele.

THE NUMBER OF WELLS has been increased from 13 to 15.

ALLELE COVERAGE:

B*78:01 to B*78:07, i.e. all the currently recognized HLA-B*78 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-July-16, release 3.1.0.

RESOLUTION IN HLA-B*78 HOMO- AND HETEROZYGOTES:

The B*78:01,78:01, B*78:01,78:02 and B*78:02,78:02 genotypes give rise to unique amplification patterns.

INFLUENCE ON THE INTERPRETATION OF HLA-B*78 SUBTYPINGS BY NON-HLA-B*78 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.