

HLA-B\*78  
101.551-06 – licensed for PCR  
101.551-06u – not licensed for PCR  
Lot No.: Y33

1

## Release Note - *Olerup* SSP™ HLA-B\*78 – Lot No. Y33

Product number:	101.551-06 – licensed for PCR 101.551-06u – <u>not</u> licensed for PCR
Lot number:	Y33
Expiry date:	2009-September-01
Number of tests:	6
Number of tubes per test:	12

### CHANGES COMPARED TO THE PREVIOUS HLA-B\*78 LOT X05):

The HLA-B\*78 primer set is unchanged.

THE NUMBER OF WELLS is unchanged.

### ALLELE COVERAGE:

B\*7801 to B\*7805, i.e. all the currently recognized HLA-B\*78 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2007-October-05, release 2.19.0.

### RESOLUTION IN HLA-B\*78 HOMO- AND HETEROZYGOTES:

The B\*7801,7801, B\*7801,7802 and B\*7802,7802 genotypes give rise to unique amplification patterns.

### INFLUENCE ON THE INTERPRETATION OF HLA-B\*78 SUBTYPINGS BY NON-HLA-B\*78 ALLELES:

None frequently occurring.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.