

101.553-06 – including Taq polymerase

101.553-06u – without Taq polymerase

Lot No.: 18E

**Olerup SSP<sup>®</sup> HLA-B\*81**

Product number: 101.553-06 – including Taq polymerase  
101.553-06u – without Taq polymerase  
Lot number: 18E  
Expiry date: 2010-March-01  
Number of tests: 6  
Number of wells per test: 5

**CHANGES COMPARED TO THE PREVIOUS HLA-B\*81 LOT (X11):**

Well	5'-primer	3'-primer	rationale
5	New	New	New primer pair for the B*8103 allele.

THE NUMBER OF WELLS has been increased from 4 to 5.

**ALLELE COVERAGE:**

B\*8101 to B\*8103, i.e. all the currently recognized HLA-B\*81 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2008-January-11, release 2.20.0.

**RESOLUTION IN HLA-B\*81 HOMO- AND HETEROZYGOTES:**

Excellent.

**INFLUENCE ON THE INTERPRETATION OF HLA-B\*81 SUBTYPINGS BY NON-HLA-B\*81 ALLELES:**

None frequently occurring.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.