

## **Olerup SSP<sup>®</sup> HLA-B\*81**

**Product number:** 101.553-06 – including Taq polymerase  
101.553-06u – without Taq polymerase  
**Lot number:** 51G  
**Expiry date:** 2011-November-01  
**Number of tests:** 6  
**Number of wells per test:** 6

### **CHANGES COMPARED TO THE PREVIOUS HLA-B\*81 LOT (56F):**

The HLA-B\*81 primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

B\*8101 to B\*8104N, i.e. all the currently recognized HLA-B\*81 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2009-October-19, release 2.27.0.

### **RESOLUTION IN HLA-B\*81 HOMO- AND HETEROZYGOTES:**

Excellent.

### **INFLUENCE ON THE INTERPRETATION OF HLA-B\*81 SUBTYPINGS BY NON-HLA-B\*81 ALLELES:**

None frequently occurring.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No suggestions received.