

101.553-06 – including Taq polymerase

101.553-06u – without Taq polymerase

Lot No.: 99K

Olerup SSP® HLA-B*81

Product number: 101.553-06 – including Taq polymerase
101.553-06u – without Taq polymerase
Lot number: 99K
Expiry date: 2013-August-01
Number of tests: 6
Number of wells per test: 7

CHANGES COMPARED TO THE PREVIOUS HLA-B*81 LOT (51G):

The primers of the wells detailed below have been exchanged, added or modified compared to the previous lot.

| Well | 5'-primer | 3'-primer | rationale |
|------|-----------|-----------|---|
| 7 | New | New | New primer pair for the B*81:05 allele. |

THE NUMBER OF WELLS is increased from 6 to 7.

ALLELE COVERAGE:

B*81:01 to B*81:05, i.e. all the currently recognized HLA-B*81 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.

RESOLUTION IN HLA-B*81 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*81 SUBTYPINGS BY NON-HLA-B*81 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.