

Olerup SSP® HLA-B*82

Product number: 101.552-06 – including *Taq* pol.
101.552-06u – without *Taq* pol.
Lot number: 27K
Expiry date: 2012-July-01
Number of tests: 6
Number of wells per test: 5

CHANGES COMPARED TO THE PREVIOUS HLA-B*82 LOT (57F):

Well	5'-primer	3'-primer	rationale
2	-	-	Exchanged positive control primer pair.
5	New	New	New primer pair for the B*82:03 allele.

THE NUMBER OF WELLS is increased from 4 to 5.

ALLELE COVERAGE:

B*82:01 to B*82:03, i.e. all the currently recognized HLA-B*82 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-April-01, release 3.0.0.

RESOLUTION IN HLA-B*82 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*82 SUBTYPINGS BY NON-HLA-B*82 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.