

Olerup SSP[®] HLA-B*82

Product number:	101.552-06 – including <i>Taq</i> pol. 101.552-06u – without <i>Taq</i> pol.
Lot number:	33M
Expiry date:	2013-November-01
Number of tests:	6
Number of wells per test:	5

CHANGES COMPARED TO THE PREVIOUS HLA-B*82 LOT (27K):

The HLA-B*82 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*82:01 to B*82:03, i.e. all the currently recognized HLA-B*82 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2011-April-08, release 3.4.0.

RESOLUTION IN HLA-B*82 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*82 SUBTYPINGS BY NON-HLA-B*82 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.