

Olerup SSP[®] HLA-B*82

Product number: 101.552-06 – including *Taq* pol.
101.552-06u – without *Taq* pol.
Lot number: 58R
Expiry date: 2015-April-01
Number of tests: 6
Number of wells per test: 5

CHANGES COMPARED TO THE PREVIOUS HLA-B*82 LOT (33M):

The HLA-B*82 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*82:01 to B*82:03, i.e. all the currently recognized HLA-B*82 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2012-July-12, release 3.9.0.

RESOLUTION IN HLA-B*82 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*82 SUBTYPINGS BY NON-HLA-B*82 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.