

## Olerup SSP® HLA-B\*82

**Product number:** 101.552-06 – including *Taq* pol.  
101.552-06u – without *Taq* pol.  
**Lot number:** 3G6  
**Expiry date:** 2021-01-01  
**Number of tests:** 6  
**Number of wells per test:** 5+1

### CHANGES COMPARED TO THE PREVIOUS HLA-B\*82 LOT (2E6):

Well	5'-primer	3'-primer	rationale
1	Modified	-	5'-primer modified for improved HLA-specific amplification.

THE NUMBER OF WELLS is unchanged.

### ALLELE COVERAGE:

B\*82:01 to B\*82:03, i.e. all the currently recognized HLA-B\*82 alleles, give rise to unique amplification patterns; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2018-April-16, release 3.32.0<sup>1</sup>.

<sup>1</sup>Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

### RESOLUTION IN HLA-B\*82 HOMO- AND HETEROZYGOTES:

Excellent.

### INFLUENCE ON THE INTERPRETATION OF HLA-B\*82 SUBTYPINGS BY NON-HLA-B\*82

#### ALLELES:

None of importance.

### MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

Primer mix 1 has been modified for improved HLA-specific amplification.