

HLA-B*82
101.552-06 – licensed for PCR
101.552-06u – not licensed for PCR
Lot No.: **Y53**

1

Release Note - *Olerup* SSP™ HLA-B*82 – Lot No. **Y53**

Product number: 101.552-06 – licensed for PCR
101.552-06u – not licensed for PCR
Lot number: Y53
Expiry date: 2009-October-01
Number of tests: 6
Number of tubes per test: 4

CHANGES COMPARED TO THE PREVIOUS HLA-B*82 LOT (V98):

Tube	5'-primer	3'-primer	rationale
3	Exchanged	Modified	Increased resolution.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

B*8201 to B*8212, i.e. all the currently recognized HLA-B*82 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2007-October-05, release 2.19.0.

RESOLUTION IN HLA-B*82 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-B*82 SUBTYPINGS BY NON-HLA-B*82 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.