

Release Note - Olerup SSP™ HLA-C low resolution – Lot No. Y22

Product number: 101.601-24/12 – licensed for PCR
101.601-24u/12u – not licensed for PCR
Lot number: Y22
Expiry date: 2009-August-01
Number of tests: 24 tests – Product No. 101.601-24
12 tests – Product No. 101.601-12
Number of tubes per test: 23 + 1

CHANGES COMPARED TO THE PREVIOUS HLA-Cw LOW RESOLUTION LOT (V97):

Tube	5'-primer	3'-primer	rationale
3	-	-	Changed to the 800 bp internal positive control primer pair to decrease primer dimer formation tendency.
5	Added	Added	Primer added for the Cw*1214 allele.
8	-	Added	Primer added for the Cw*0603 allele.
11	Modified	-	Increased yield of specific PCR product.
16	Added	Added	Primer added for the Cw*0511 allele.
17	Added	-	Primer added for the Cw*1520 allele.
18	-	Added	Primer added for the Cw*0427 allele.
22	Added	Added	Primer added for the Cw*0511 allele.
23	Added	-	Primer added for the Cw*1610 allele.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

All the HLA-Cw alleles, recognized by the HLA Nomenclature Committee in October 2007 will be amplified by the primers in the HLA-Cw low resolution SSP kit. The HLA-Cw alleles will be grouped according to the first two digits of the allele designations; www.ebi.ac.uk/imgt/hla, 2007-October-05, release 2.19.0.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.