

## **Olerup SSP<sup>®</sup> HLA-C\*04 Add-on**

<b>Product number:</b>	<b>101.862-12 – including <i>Taq</i> polymerase 101.862-12u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>56R</b>
<b>Expiry date:</b>	<b>2015-April-01</b>
<b>Number of tests:</b>	<b>12</b>
<b>Number of wells per test:</b>	<b>4</b>

**THE NUMBER OF WELLS** is 4 in an 8 well plate.

### **ALLELE COVERAGE:**

The HLA-C\*04:01, 04:09N, 04:30 and 04:82 alleles give different patterns in the HLA-C\*04 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2012-July-12, release 3.9.0.

<sup>1</sup>The HLA-C\*04 add-on kit cannot distinguish the silent mutations in the C\*04:01:01:01-04:01:41 alleles.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Good.

### **INFLUENCE ON THE INTERPRETATION OF HLA-C\*04 ADD-ON SUBTYPINGS BY NON-HLA-C\*04 ALLELES:**

None frequently occurring.