

Olerup SSP[®] HLA-C*04 Add-on

Product number:	101.862-12 – including <i>Taq</i> polymerase 101.862-12u – without <i>Taq</i> polymerase
Lot number:	49S
Expiry date:	2015-December-01
Number of tests:	12
Number of wells per test:	4

THE NUMBER OF WELLS is 4 in an 8 well plate.

ALLELE COVERAGE:

The HLA-C*04:01, 04:09N, 04:30 and 04:82 alleles give different patterns in the HLA-C*04 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2013-April-17, release 3.12.0.

¹The HLA-C*04 add-on kit cannot distinguish the silent mutations in the C*04:01:01:01-04:01:51 alleles.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*04 ADD-ON SUBTYPINGS BY NON-HLA-C*04 ALLELES:

None frequently occurring.