

101.862-12 – including *Taq* polymerase101.862-12u – without *Taq* polymerase

Lot No.: 53Y

Olerup SSP® HLA-C*04 Add-on

Product number: 101.862-12 – including *Taq* polymerase
101.862-12u – without *Taq* polymerase

Lot number: 53Y

Expiry date: 2017-November-01

Number of tests: 12

Number of wells per test: 4+1

CHANGES COMPARED TO THE PREVIOUS HLA-C*04 ADD-ON LOT (92X):

The HLA-C*04 Add-on primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The HLA-C*04:01, 04:09N, 04:30 and 04:82 alleles give different patterns in the HLA-C*04 Add-on kit^{1,2}; www.ebi.ac.uk/imgt/hla, 2015-January-19, release 3.19.0.

¹The HLA-C*04 add-on kit cannot distinguish the silent mutations in the C*04:01:01:01-04:01:66 alleles.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*04 ADD-ON SUBTYPINGS BY NON-HLA-C*04 ALLELES:

None frequently occurring.