

101.862-12 – including *Taq* polymerase101.862-12u – without *Taq* polymeraseLot No.: **6E1****Olerup SSP® HLA-C*04 Add-on**

Product number:	101.862-12 – including <i>Taq</i> polymerase 101.862-12u – without <i>Taq</i> polymerase
Lot number:	6E1
Expiry date:	2019-05-01
Number of tests:	12
Number of wells per test:	4+1

The HLA-C*04 Add-on primer set is unchanged compared to the previous *Olerup SSP® HLA-C*04 Add-on (Lot No. 8D9)*.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The HLA-C*04:01, 04:09N, 04:30 and 04:82 alleles give different patterns in the HLA-C*04 Add-on kit^{1,2}; www.ebi.ac.uk/imgt/hla, 2016-July-14, release 3.25.0.

¹The HLA-C*04 add-on kit cannot distinguish the silent mutations in the C*04:01:01:01-04:01:75 alleles.

²Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*04 ADD-ON SUBTYPINGS BY NON-HLA-C*04 ALLELES:

None frequently occurring.