

HLA-C\*04:09N  
101.861-12 – including *Taq* polymerase  
101.861-12u – without *Taq* polymerase  
Lot No.: **52S**

Release Note

Page 1 of 1

## **Olerup SSP<sup>®</sup> HLA-C\*04:09N**

<b>Product number:</b>	<b>101.861-12 – including <i>Taq</i> polymerase 101.861-12u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>52S</b>
<b>Expiry date:</b>	<b>2015-November-01</b>
<b>Number of tests:</b>	<b>12</b>
<b>Number of wells per test:</b>	<b>2</b>

### **CHANGES COMPARED TO THE PREVIOUS HLA-C\*04:09N LOT (41R):**

The HLA-C\*04:09N primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is unchanged.

### **ALLELE COVERAGE:**

The C\*04:09N allele gives rise to a unique amplification pattern;  
[www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2013-January-11, release 3.11.0.

### **RESOLUTION IN HLA-C\*04:09N HOMO- AND HETEROZYGOTES:**

Good.

### **INFLUENCE ON THE INTERPRETATION OF HLA-C\*04:09N SUBTYPINGS BY NON-HLA-C\*04 ALLELES:**

None of importance.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.