

101.861-12 – including *Taq* polymerase101.861-12u – without *Taq* polymerase

Lot No.: 51X

**Olerup SSP® HLA-C\*04:09N**

**Product number:** 101.861-12 – including *Taq* polymerase  
101.861-12u – without *Taq* polymerase

**Lot number:** 51X

**Expiry date:** 2017-May-01

**Number of tests:** 12

**Number of wells per test:** 2+1

**CHANGES COMPARED TO THE PREVIOUS HLA-C\*04:09N LOT (65V):**

Well	5'-primer	3'-primer	rationale
3	-	-	Updated negative control.

**THE NUMBER OF WELLS** is unchanged.

**ALLELE COVERAGE:**

The C\*04:09N allele gives rise to a unique amplification pattern; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2014-October-10, release 3.18.0.

**RESOLUTION IN HLA-C\*04:09N HOMO- AND HETEROZYGOTES:**

Good.

**INFLUENCE ON THE INTERPRETATION OF HLA-C\*04:09N SUBTYPINGS BY NON-HLA-C\*04 ALLELES:**

None of importance.

**MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.