

101.861-12 – including *Taq* polymerase

101.861-12u – without *Taq* polymerase

Lot No.: **65V**

## **Olerup SSP® HLA-C\*04:09N**

**Product number:** 101.861-12 – including *Taq* polymerase  
101.861-12u – without *Taq* polymerase

**Lot number:** 65V

**Expiry date:** 2016-October-01

**Number of tests:** 12

**Number of wells per test:** 2+1

### **CHANGES COMPARED TO THE PREVIOUS HLA-C\*04:09N LOT (52S):**

Well	5'-primer	3'-primer	rationale
3	-	-	Negative Control

**THE NUMBER OF WELLS** is increased from 2 to 3 wells.

### **ALLELE COVERAGE:**

The C\*04:09N allele gives rise to a unique amplification pattern;  
[www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2014-January-17, release 3.15.0.

### **RESOLUTION IN HLA-C\*04:09N HOMO- AND HETEROZYGOTES:**

Good.

### **INFLUENCE ON THE INTERPRETATION OF HLA-C\*04:09N SUBTYPINGS BY NON-HLA-C\*04 ALLELES:**

None of importance.

### **MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:**

No comments received.