

101.861-12 – including *Taq* polymerase

101.861-12u – without *Taq* polymerase

Lot No.: **5D4**

Olerup SSP® HLA-C*04:09N

Product number:	101.861-12 – including <i>Taq</i> polymerase 101.861-12u – without <i>Taq</i> polymerase
Lot number:	5D4
Expiry date:	2018-07-01
Number of tests:	12
Number of wells per test:	2+1

CHANGES COMPARED TO THE PREVIOUS HLA-C*04:09N LOT (85Y):

Well	5'-primer	3'-primer	rationale
3	-	-	Updated negative control.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The C*04:09N allele gives rise to a unique amplification pattern; www.ebi.ac.uk/imgt/hla, 2015-October-10, release 3.22.1.

RESOLUTION IN HLA-C*04:09N HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*04:09N SUBTYPINGS BY NON-HLA-C*04 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

None.