

HLA-C*04:09N
101.861-12 – including *Taq* polymerase
101.861-12u – without *Taq* polymerase
Lot No.: **5F7**

Olerup SSP[®] HLA-C*04:09N

Product number:	101.861-12 – including <i>Taq</i> polymerase 101.861-12u – without <i>Taq</i> polymerase
Lot number:	5F7
Expiry date:	2021-09-01
Number of tests:	12
Number of wells per test:	2+1

The HLA-C*04:09N primer set is unchanged compared to the previous *Olerup* SSP[®] HLA-C*04:09N (Lot No. 3E3).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The C*04:09N allele gives rise to a unique amplification pattern; www.ebi.ac.uk/imgt/hla, 2017-April-13, release 3.28.0.

RESOLUTION IN HLA-C*04:09N HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*04:09N SUBTYPINGS BY NON-HLA-C*04 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

None.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to extension of shelf-life.