

101.861-12 – including *Taq* polymerase

101.861-12u – without *Taq* polymerase

Lot No.: **4L4**

Olerup SSP® HLA-C*04:09N

Product number:	101.861-12 – including <i>Taq</i> polymerase 101.861-12u – without <i>Taq</i> polymerase
Lot number:	4L4
Expiry date:	2024-08-01
Number of tests:	12
Number of wells per test:	2+1

CHANGES COMPARED TO THE PREVIOUS OLERUP SSP® HLA-C*04:09N

LOT (5F7):

The HLA-C*04:09N primer set is unchanged compared to the previous *Olerup* SSP® HLA-C*04:09N (Lot No. 5F7).

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The C*04:09N allele gives rise to a unique amplification pattern; www.ebi.ac.uk/imgt/hla, 2020-April-2020, release 3.40.0.

RESOLUTION IN HLA-C*04:09N HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*04:09N SUBTYPINGS BY NON-HLA-C*04 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

None.