

Olerup SSP[®] HLA-C*04:09N

Product number:	101.861-12 – including <i>Taq</i> polymerase 101.861-12u – without <i>Taq</i> polymerase
Lot number:	41R
Expiry date:	2015-March-01
Number of tests:	12
Number of wells per test:	2

CHANGES COMPARED TO THE PREVIOUS HLA-C*04:09N LOT (15N):

The HLA-C*04:09N primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The C*04:09N allele gives rise to a unique amplification pattern;
www.ebi.ac.uk/imgt/hla, 2012-July-12, release 3.9.0.

RESOLUTION IN HLA-C*04:09N HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*04:09N SUBTYPINGS BY NON-HLA-C*04 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.