

Olerup SSP® HLA-C*07 Add-on

Product number:	101.863-12 – including <i>Taq</i> polymerase 101.863-12u – without <i>Taq</i> polymerase
Lot number:	44S
Expiry date:	2015-December-01
Number of tests:	12
Number of wells per test:	5

CHANGES COMPARED TO THE PREVIOUS HLA-C*07 ADD-ON LOT (26N):

The HLA-C*07 Add-on primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The HLA-C*07:06 and C*07:18 and the C*07:01, the C*07:50 and the C*07:02 and the C*07:11 and the C*07:04 alleles are uniquely separated from each other in the HLA-C*07 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2013-April-17, release 3.12.0.

¹The HLA-C*07 Add-on kit cannot distinguish the silent mutations in the C*07:01:01-07:01:33 alleles, the 07:02:01:01-07:02:41 alleles or the C*07:04:01-C*07:04:08 alleles.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*07 ADD-ON SUBTYPINGS BY NON-HLA-C*07 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.