

Olerup SSP® HLA-C*08 Add-on

Product number:	101.862-12 – including <i>Taq</i> polymerase 101.862-12u – without <i>Taq</i> polymerase
Lot number:	27N
Expiry date:	2014-August-01
Number of tests:	12
Number of wells per test:	2

THE NUMBER OF WELLS is 2 in an 8 well plate.

ALLELE COVERAGE:

The HLA-C*08:01, and 08:22 alleles give different amplification patterns in the HLA-C*08 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2012-January-12, release 3.7.0.

¹The HLA-C*08 Add-on kit cannot distinguish the silent mutations in the C*08:01:01-08:01:05 alleles.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*08 ADD-ON SUBTYPINGS BY NON-HLA-C*08 ALLELES:

None frequently occurring.