

## Olerup SSP<sup>®</sup> HLA-C\*08 Add-on

**Product number:** 101.864-12 – including *Taq* polymerase  
101.864-12u – without *Taq* polymerase  
**Lot number:** 45S  
**Expiry date:** 2015-December-01  
**Number of tests:** 12  
**Number of wells per test:** 2

THE NUMBER OF WELLS is 2 in an 8 well plate.

### ALLELE COVERAGE:

The HLA-C\*08:01, and 08:22 alleles give different amplification patterns in the HLA-C\*08 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2013-April-17, release 3.12.0.

<sup>1</sup>The HLA-C\*08 Add-on kit cannot distinguish the silent mutations in the C\*08:01:01-08:01:09 alleles.

### RESOLUTION IN HOMO- AND HETEROZYGOTES:

Good.

### INFLUENCE ON THE INTERPRETATION OF HLA-C\*08 ADD-ON SUBTYPINGS BY NON-HLA-C\*08 ALLELES:

None frequently occurring.