

101.864-12 – including *Taq* polymerase**101.864-12u – without *Taq* polymerase****Lot No.: 0H9**

Olerup SSP[®] HLA-C*08 Add-on

Product number: 101.864-12 – including *Taq* polymerase
101.864-12u – without *Taq* polymerase

Lot number: 0H9

Expiry date: 2022-10-01

Number of tests: 12

Number of wells per test: 4+1

CHANGES COMPARED TO THE PREVIOUS HLA-C*08 ADD-ON LOT (6E7):

The HLA-C*08 Add-on primer set is unchanged compared to the previous *Olerup SSP[®] HLA-C*08 Add-on (Lot No. 6E7)*.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

The HLA-C*08:01, 08:22, 08:99 and 08:102 alleles give different amplification patterns in the HLA-C*08 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2018-July-11, release 3.33.0.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*08 ADD-ON SUBTYPINGS BY NON-HLA-C*08 ALLELES:

None frequently occurring.

Changes in revision R01 compared to R00:

1. The expiration date has been altered due to shelf-life extension.