

## **Olerup SSP<sup>®</sup> HLA-C\*15 Add-on**

<b>Product number:</b>	<b>101.865-12 – including <i>Taq</i> polymerase</b> <b>101.865-12u – without <i>Taq</i> polymerase</b>
<b>Lot number:</b>	<b>47S</b>
<b>Expiry date:</b>	<b>2015-December-01</b>
<b>Number of tests:</b>	<b>12</b>
<b>Number of wells per test:</b>	<b>2</b>

### **CHANGES COMPARED TO THE PREVIOUS HLA-C\*15 ADD-ON LOT (28N):**

The HLA-C\*15 Add-on primer set is unchanged compared to the previous lot.

**THE NUMBER OF WELLS** is 2 in an 8 well plate.

### **ALLELE COVERAGE:**

The HLA-C\*15:05 and 15:29 alleles give different patterns in the HLA-C\*15 Add-on kit<sup>1</sup>; [www.ebi.ac.uk/imgt/hla](http://www.ebi.ac.uk/imgt/hla), 2013-April-17, release 3.12.0.

<sup>1</sup>The HLA-C\*15 Add-on kit cannot distinguish the silent mutations in the C\*15:05:01-15:05:08 alleles.

### **RESOLUTION IN HOMO- AND HETEROZYGOTES:**

Good.

### **INFLUENCE ON THE INTERPRETATION OF HLA-C\*15 ADD-ON SUBTYPINGS BY NON-HLA-C\*15 ALLELES:**

None.