

Olerup SSP[®] HLA-C*15 Add-on

Product number:	101.865-12 – including <i>Taq</i> polymerase 101.865-12u – without <i>Taq</i> polymerase
Lot number:	28N
Expiry date:	2014-August-01
Number of tests:	12
Number of wells per test:	2

THE NUMBER OF WELLS is 2 in an 8 well plate.

ALLELE COVERAGE:

The HLA-C*15:05 and 15:29 alleles give different patterns in the HLA-C*15 Add-on kit¹; www.ebi.ac.uk/imgt/hla, 2012-January-12, release 3.7.0.

¹The HLA-C*15 Add-on kit cannot distinguish the silent mutations in the C*15:05:01-15:05:06 alleles.

RESOLUTION IN HOMO- AND HETEROZYGOTES:

Good.

INFLUENCE ON THE INTERPRETATION OF HLA-C*15 ADD-ON SUBTYPINGS BY NON-HLA-C*15 ALLELES:

None.