

101.628-06 – including *Taq* polymerase

101.628-06u – without *Taq* polymerase

Lot No.: **4H9**

Olerup SSP[®] HLA-C*17

Product number:	101.628-06 – including <i>Taq</i> polymerase 101.628-06u – without <i>Taq</i> polymerase
Lot number:	4H9
Expiry date:	2021-09-01
Number of tests:	6
Number of wells per test:	15+1

CHANGES COMPARED TO THE PREVIOUS HLA-C*17 LOT (2G9):

The HLA-C*17 primer set is unchanged compared to the previous lot.

THE NUMBER OF WELLS is unchanged.

ALLELE COVERAGE:

C*17:01 to C*17:41, i.e. all the currently recognized HLA-C*17 alleles, will be amplified by the primers in the HLA-C*17 subtyping kit; www.ebi.ac.uk/imgt/hla, 2019-January-23¹, release 3.35.0.

The HLA-C*17 kit enables separation of the confirmed HLA-C*17 alleles as listed in the IMGT/HLA database 3.24.0. An HLA allele is listed as confirmed by IMGT/HLA if it has been sequenced by more than a single laboratory or from multiple sources.

The HLA-C*17 kit also enables identification of null and alternatively expressed alleles.

¹Alleles that have been deleted from or renamed in the official WHO HLA Nomenclature up to and including the last IMGT/HLA database release can be retrieved from web page <http://hla.alleles.org/alleles/deleted.html>.

RESOLUTION IN HLA-C*17 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-C*17 SUBTYPINGS BY NON-HLA-C*17 ALLELES:

None of importance.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No comments received.