

Olerup SSP[®] HLA-C*18

Product number: 101.629-06 – including *Taq* polymerase
101.629-06u – without *Taq* polymerase
Lot number: 97K
Expiry date: 2013-August-01
Number of tests: 6
Number of wells per test: 6

CHANGES COMPARED TO THE PREVIOUS HLA-C*18 LOT (97K):

The primers of the wells detailed below has been exchanged, added or modified compared to the previous lot.

Well	5'-primer	3'-primer	rationale
6	New	New	New primer pair for the C*18:04 allele.

THE NUMBER OF WELLS is increased from 5 to 6.

ALLELE COVERAGE:

C*18:01 to C*18:04, i.e. all the currently recognized HLA-C*18 alleles, give rise to unique amplification patterns; www.ebi.ac.uk/imgt/hla, 2010-October-15, release 3.2.0.

RESOLUTION IN HLA-C*18 HOMO- AND HETEROZYGOTES:

Excellent.

INFLUENCE ON THE INTERPRETATION OF HLA-C*18 SUBTYPINGS BY NON-HLA-C*18 ALLELES:

None frequently occurring.

MODIFICATIONS MADE DUE TO COMMENTS FROM CUSTOMERS:

No suggestions received.